

Aesculap[®] Columbus[®] Rationale

One world – One knee



Aesculap Orthopaedics



Contents

	Page
List of key words	4
List of abbreviations	5
1. Introduction	7
2. The Columbus® Portfolio	8
2.1 Overview	8
2.1.1 PCL-retaining fixed gliding surface Columbus® CR DD	9
2.1.2 PCL-retaining rotating gliding surface Columbus® RP	9
2.1.3 PCL-sacrificing fixed gliding surface Columbus® PS	10
2.1.4 PCL-sacrificing fixed gliding surface Columbus® UC	11
2.1.5 PCL-sacrificing rotating gliding surface Columbus® UCR	12
2.1.6 Patella	14
2.1.7 Tibia augment and extension stems	14
2.1.8 Columbus® Revision	15
2.2 Implant sizes	16
3. The Columbus® Design	18
3.1 Patellofemoral kinematics	18
3.2 Tibiofemoral kinematics	20
3.3 Abrasion	21
3.4 Kinematic gait analysis	24
3.5 Stability	26
4. Selected materials	28
4.1 Surface properties	28
4.2 Plasmapore® cementless fixation	28
4.3 Beta polyethylene	29
4.4 AS coating	30
5. Versatile instruments	31
5.1 Overview	31
5.2 OrthoPilot® navigation system	31
5.3 Columbus® MIOS® – Minimally Invasive Orthopaedic Solutions	33
6. Clinical outcomes	34
7. Important parameters for Columbus® implants	38
8. Implant selection matrix Columbus® Primary	41
9. Implant selection matrix Columbus® AS	43
10. Implant selection matrix Columbus® Revision	45
Literature	48

List of key words

Key word	Page		
Abrasion	21	Materials	28 et seq.
Allergy	30	Maximum flexion angle	20
Anatomical fit	16, 17, 18, 19	Mechanical axis	31
Anterior knee pain	34	Meniscus components	9, 10, 11
Augmentation	14	Minimally invasive surgical technique	32
Beta polyethylene	29	Navigation system	31
Cement adhesion	33	OrthoPilot®	31
Cementless implants	28	Oxford Score	35
Competition systems	23	Patella	14
Cross linking	29	Plasmapore® coating	28
Defect	14, 15	PMMA nanobond coating	33
Delamination	21	Posterior slope	20
Extension stem	14, 15	Post-operative outcomes	34 et seq.
Femoral design	18, 19	Prostheses register	34
Femoral radius	18	Revision system	15
Fixation	28	Revision rate	34
Fixed platform	27	X-ray images	36
Gait cycle analysis	24, 25	Rotating platform	35
Posterior cruciate ligament	8, 9, 10, 11, 12, 13, 24	Scientific results	34 et seq.
Highly cross-linked polyethylene	29	Stability	26
Implant size	16, 22	Studies	32
Implant dimensions	38, 39	Tibia	17, 27
Indication	8	Tibial design	17, 27
Instrumentation	7, 12	Survival rate	23, 34 et seq.
Instability	10, 11, 12, 15	Varus / Valgus	31, 36
Jumping distance	10, 11		
Kinematics	20, 24, 25		
Clinical outcomes	34 et seq.		
Compatibility	8, 13, 30		
Congruence	22		
Contact surface	21		
KSS	34, 35, 36		

List of abbreviations

A	Anterior
AS	Advanced Surface
CR	Cruciate Retaining
CRA	Cruciate Retaining Augmentation
PSA	Posterior Stabilized Augmentation
DD	Deep Dish
F	Fix
HC	High Constraint
HXLPE	Highly Cross-Linked Polyethylene
PCL	Posterior Cruciate Ligament
KSS-F	Knee Society Score – Function
KSS-K	Knee Society Score – Clinical
L	Lateral
M	Medial
MC	Medium Constraint
MIOS®	Minimally Invasive Orthopedic Solutions
N	Narrow
OXF	Oxford-Score
P	Posterior
PE	Polyethylene
PS	Posterior Stabilized
RP	Rotating Platform
UC	Ultra Congruent
UCR	Ultra Congruent Rotating
UHMWPE	Ultra High Molecular Weight Polyethylene

One world – One knee



1. Introduction

Although revision rates in knee replacement procedures have been continuously decreasing over recent years, there has been little improvement in patient satisfaction, with the number of dissatisfied patients following knee arthroplasty remaining at around 14 %.¹ The main complications following knee replacement surgery include unnatural tibiofemoral kinematics, an abnormal patellar motion, polyethylene abrasion and a limited degree of flexion.^{2,3,4} The Columbus® knee prosthesis was developed to resolve exactly these problems with its advanced implant design and versatile options in the selection of implants and instruments, and to provide a good surgical outcome for both the patient and surgeon.

Superior design features of the Columbus®:

- Natural kinematics
- Excellent abrasion rate
- High bending capability
- Outstanding patellar motion



Columbus® knee prosthesis

2. The Columbus® Portfolio

Variety – Discover individuality

2.1 Overview

The Columbus® Knee System offers a suitable solution for nearly all patient needs. The choice is yours: The combinability of Columbus® allows an intraoperative change from a PCL-retaining (CR DD, RP) treatment to a PCL-sacrificing (UC, UCR, PS) treatment. The implant options vary with regard to the degree of guidance and stabilization, and they therefore enable the appropriate treatment of nearly all conditions of the ligaments. Minor defects can be treated in the primary

region with tibial plates, tibial extension stems and gliding surfaces up to 20 mm in height. For severe bone defects and unstable ligament conditions, the Columbus® Revision System provides extensive treatment options.

Discover Compatibility




Figure 1: Columbus® product overview

2.1.1 PCL-retaining fixed gliding surface Columbus® CR DD

Thanks to its high stability, the PCL-retaining, fixed gliding surface CR DD is also suitable for difficult ligament conditions, such as a ruptured PCL, and at the same time offers natural kinematics.⁵ The CR DD design is characterized by a deep concavity in the sagittal plane, enabling a high level of congruence to be achieved between the femoral components and the gliding surface. Stability is therefore increased particularly in the AP direction compared with the flatter gliding surface, without restricting bending capability.




Figure 2: Sagittal section of the Columbus® CR DD meniscus components compared with the flatter CR meniscus design.

2.1.2 PCL-retaining rotating gliding surface Columbus® RP

The Columbus® RP variant allows a physiological rotational movement with a rotational freedom of $\pm 10^\circ$. An optimum tibial coverage can also be achieved more easily. The RP design has been in successful use since 2002, offering high protection against luxation and an almost natural sequence of movements in the joint.



Figure 3: Columbus® RP gliding surface with tibial plateau.

2. The Columbus® Portfolio

Variety – Discover individuality

2.1.3 PCL-sacrificing fixed gliding surface Columbus® PS

In the event of PCL insufficiency, the Columbus® PS prosthesis takes on the role of the PCL and assures a dorsal rollback through its post-cam mechanism. With the Columbus® PS knee prosthesis, the "Jumping Distance" (JP) continuously increases with increasing flexion, enabling the increased risk of luxation with large flexion angles to be avoided.⁶

The limited translational and rotational possibilities increase with the size of the femur, allowing freedom of movement to be adapted to the anatomical conditions.



JP	30°	60°	90°	120°
F1	13.7	15.2	15.1	15.2
F2	12.8	14.4	14.4	14.4
F3	11.9	13.5	13.7	13.5
F4	11.2	13	13.4	13.1
F5	9.5	12.3	12.9	12.1
F6	8.6	11.7	12.7	11.9
F7	7.9	11.2	12.4	11.7

Figure 4: Jumping Distance of the Columbus® PS prosthesis in mm.


Columbus® PS	Translation AP in mm	Translation ML in mm	Rotation Int. / Ext.
F1	5.8 mm	± 1 mm	± 10°
F2	5.9 mm	± 1.25 mm	± 10°
F3	6 mm	± 1.75 mm	± 12°
F4	6.1 mm	± 1.75 mm	± 12°
F5	6.2 mm	± 2 mm	± 14°
F6	6.2 mm	± 2.5 mm	± 16°
F7	6.3 mm	± 3.5 mm	± 20°
F8	6.5 mm	± 3.5 mm	± 20°

Translational and rotational freedoms of the Columbus® PS prosthesis.

2.1.4 PCL-sacrificing fixed gliding surface Columbus® UC

Treatment with the PCL-sacrificing Columbus® UC is easier on the bone than the PS variant, as no femoral box needs to be prepared for a post-cam mechanism.

The ventrally elevated edge also offers high stability in the AP direction when resecting the PCL.⁷



Benefits of the Columbus® UC:

- Additional option for cruciate ligament treatment.
- Bone sparing, as no bone needs to be removed for the PS box.
- Time saving, as there is no PS box preparation.
- High AP stability due to the ventrally elevated edge.
- Low risk of patellar impingement due to the ventrally flattened meniscus design.
- The decision regarding which Columbus® UC variant to use does not have to be made until after the final implantation of the metal components.

Jumping Distance

	T1	T2	T3	T4	T5
A	7.6	8	8.7	9.2	10.2
B	11.7	11.6	12.1	12	12.8

Figure 5: Jumping Distance of the Columbus® UC gliding surface compared with the Columbus® DD gliding surface in mm.



Columbus® UC gliding surface


2. The Columbus® Portfolio

Variety – Discover individuality

2.1.5 PCL-sacrificing rotating gliding surface Columbus® UCR


The Columbus® UCR knee prosthesis enables a rotating gliding surface to be used for PCL resections:

- Rotational freedom of $\pm 20^\circ$ provides great tibial coverage.
- The tibial overhang does not cause any adverse impingement due to the limited rotational freedom and a bevel of the gliding surface in the dorsal region (Fig. 6).
- The tibial implant stem length has been standardized at 28 mm for all tibia sizes.
- The rounded rotation stop ensures high stability.
- The ventral elevation of the PE components leads to high AP stability (Fig. 7).
- Elevated contact surfaces reduce the abrasion rate to under 1 mg per million cycles.
- An additional instrument set provides easy instrumentation.



Size	20° Rotation		
	A (mm)	B (mm)	Surface (mm ²)
T5 (Worst Case)	3.6	0.8	75

Figure 6: Tibial overhang of the Columbus® UCR prosthesis in the worst case.



Size	T0	T1	T2	T3	T4	T5
B	18.8 mm	19.3 mm	19.7 mm	20.1 mm	20.5 mm	20.9 mm

Figure 7: Height of the ventral edge of a Columbus® UCR gliding surface with a height of 10 mm.

As all sizes are compatible for PCL-retaining treatments, the respective appropriate sizes of the prostheses can be used for both the femur and tibia.

For PCL-sacrificing prostheses, combinations of sizes are possible, but are more limited.

- Recommended: minimum abrasion, highest stability, best kinematics
- Possible: slightly increased abrasion, high stability and good kinematics
- To be avoided: increased abrasion, low stability

Compatibility
FEMUR PS / PE PS
FEMUR CR / PE UC
FEMUR CR / PE UCR

Sizes	F1	F2 F2N	F3 F3N	F4 F4N	F5 F5N	F6 F6N	F7 F7N	F8 F8N
T0/T0+								
T1/T1+								
T2/T2+								
T3/T3+								
T4/T4+								
T5								

Compatibility
FEMUR Revision / PE HC or MC

Sizes	F1	F2	F3	F4	F5	F6	F7
T0/T0+							
T1/T1+							
T2/T2+							
T3/T3+							
T4/T4+							
T5							

Figure 8: Compatibility of the PCL-sacrificing Columbus® prostheses.

2. The Columbus® Portfolio

Variety – Discover individuality

2.1.6 Patella

There are four different sizes available for replacement for the posterior surface of the patella, which can be combined with any femur size. The radius of the patella corresponds with the radius of the trochlea for all sizes.

Therefore, despite the possibility of combining all sizes with each other, there is always an optimum contact surface between the trochlea and patella. The thickness of the replacement for the posterior surface of the patella increases by 1 mm for each increasing size. The three fixation pins ensure a rotationally stable fixation and their conical form is also more sparing on the bone than a central pin.



Columbus® Patella

2.1.7 Tibia augment and extension stems

For the treatment of small tibial defects, augment at a height of 4 and 8 mm are also available as part of primary treatment. Additional stability of the tibial implant can also be achieved using the extension stems.



Columbus® augments and extension stems

2.1.8 Columbus® Revision

Severe bone defects can be compensated with the Columbus® Revision prosthesis with wedges up to 15 mm. For unstable ligament conditions, the system offers gliding surfaces in a "Medium Constraint" and "High Constraint" design to support the implant. The limited ML translation and increased varus / valgus rotational stability provides good bending capability with a lower risk of luxation. The tibial and femoral offset stems, which are available in various lengths and diameters as cementless 5° / 7° and cemented 6°, provide appropriate support for the patient in the medullary cavity.




Columbus® Revision




Tibial offset possibilities ML +/- 6 mm.




Femoral offset possibilities ML +/- 2 mm and AP up to +/- 4 mm.



4 Longitudinal grooves to prevent the risk of embolism



Cemented stem



10 Longitudinal grooves for a secure fixation in the medullary cavity (Wagner profile)



Cementless stem

Figure 9: Stem options of the Columbus® Revision System.

2. The Columbus® Portfolio

Variety – Discover individuality

2.2 Implant sizes

A comparison of the Columbus® femoral components with the femur sizes measured in the study by Yoshioka et al. indicates virtually complete coverage of all anatomical data.⁸ With 13 femur sizes (8 standard sizes and 5 narrow sizes), the Columbus® Knee System provides good coverage for both the female and male morphotypes.⁹

In ML, the narrow sizes "N" are one size smaller than the standard sizes.




Figure 10: AP / ML ratio shown according to female and male morphotypes.



	ML dimensions						
AP dimensions	F1	F2	F3	F4	F5	F6	F7
F1	F1						
F2	F2N	F2					
F3		F3N	F3				
F4			F4N	F4			
F5				F5N	F5		
F6					F6N	F6	
F7							F7
F8							F8

Figure 11: ML dimensions of Columbus® Femur Standard and in narrower "N" sizes.

Good tibial coverage is achieved by selecting a total of 9 symmetrical tibia sizes (5 standard sizes and 4 plus sizes).⁹





Figure 12: Tibial AP / ML ratio of 337 knees versus tibial AP / ML ratio of various implant systems.⁹

In AP the tibial plus sizes are 3-4 mm bigger than the standard sizes.



3. The Columbus® Design

Mobility – Discover stability

3.1 Patellofemoral kinematics

The multiple femoral radii imitate the anatomy of the human knee and ensure homogeneous extension and flexion. Tension peaks are avoided both in the patella and in the gliding surface.^{12, 13}

The Columbus® femur prosthesis has the following design features:

- The anterior 3° tapering of the femoral shield reduces the risk of anterior notching and increases the contact pressure of the implant on the bone.
- The patella radius (Rz) enables a constant patellar pressure up to 110° flexion and results in an anatomical patellar motion.^{13, 14}
- The transition radius (Ry) creates a round sagittal femoral design, helping to avoid tension in the retinaculum and patellar ligament.
- The small dorsal femoral radius (Rx) and the short posterior condyles enable high flexion up to 140°.
- Due to the rounder Columbus® femoral design, an equal collateral ligament tension is achieved in flexion and extension. Tension peaks in the PE can therefore be avoided.




Figure 13: Anatomy and geometry of the distal femur.¹³




Figure 14: Sagittal profile of the Columbus® femur prosthesis.




Figure 15: Columbus® femoral design compared with a more angular femoral design.

- The extended trochlear motion enables a surface contact between the femur and patella up to 110° flexion and therefore reduces stresses in the PE caused by point contacts.



Figure 16: Surface contact between the patella and femur up to a flexion angle of 110°.

- The repositioned trochlea helps prevent too much patellar pressure.




Figure 17: Repositioning of the patella.

- The elevated antero-lateral femoral design helps prevent subluxation of the patella, including in the early flexion stage.¹⁵




Figure 18: Elevated antero-lateral femoral design.

- According to Eckhoff, the trochlea with an anatomical valgus of 7° permits optimum patellar motion and load relief for the knee extensors.¹⁶



Figure 19: Anatomical valgus of 7° in the Columbus® trochlea.

3. The Columbus® Design

Mobility – Discover stability


3.2 Tibiofemoral kinematics

The posterior slope of 3° built into all gliding surfaces facilitates deep flexion and prevents the increase of anterior micro-movements. Bai shows that this helps further reduce tension in the PE.¹⁷

With the anatomical slope in the gliding surface, the tibial resection can be carried out with 0° slope. This considerably reduces the shearing forces between the implant and bone.



Figure 20: Posterior slope of 3° in the Columbus® gliding surface.



Columbus® Design	RP	CR DD	UC	PS	Revision MC / HC
Theoretical maximum overextension	10°	10°	10°	5°	4°
Theoretical maximum flexion	140°	140°	140°	135°	130°

Figure 21: Maximum flexion angle of the Columbus® prosthesis.

3.3 Abrasion

The highly congruent design of Columbus® increases the contact surface between the femur and PE, helping to avoid peak stresses in the polyethylene, which are responsible for delamination and early implant failure.^{20, 21, 22}





Figure 22: Tension and delamination risk in flat PE designs compared with the Columbus® DD PE Design.²⁶



Note: All knee systems marked with * are not fully comparable, as the tests were carried out with a lesser force.

Figure 23: Contact surface of various knee implant systems.^{18, 19}

3. The Columbus® Design

Mobility – Discover stability

Bartel et al show that there should be high congruence between the femur and PE above all in the medio-lateral direction to prevent high contact tension.²³

The radii of all Columbus® femur and meniscus components are identical in the frontal plane for all implant sizes to maintain congruence, even when femur and PE components of different sizes are combined with each other.



Figure 24: Optimized congruence between the femur and PE shown in the frontal plane.




Figure 25: Small femoral components combined with larger tibial components shown in the frontal plane.

The reduced abrasion rate of the Columbus® knee prostheses, of up to less than 1 mg / million cycles, reduces the risk of osteolysis, and therefore the implant loosening rate.^{24, 25}




Figure 26: Abrasion rate of the Columbus® knee implant systems according to ISO Standard 14243.^{26 - 31}

Note:

All results marked with * were carried out with ISO 14243.

3. The Columbus® Design

Mobility – Discover stability

3.4 Kinematic gait analysis

The gait cycle analysis according to ISO 14243-1 shows that, despite being highly stable, Columbus® enables natural kinematics to high extent. The AP force and internal / external torque represent the natural force effect during a gait cycle in a healthy knee. It has been established for all Columbus® implant variants that the AP force during the stance phase (0 to 0.5 seconds) barely leads to an AP shift. The high congruence of the prostheses at low flexion angles results in the AP force being largely absorbed by the prostheses, in addition to a reduction of contact tension.

In the swing phase (0.5 to 1 seconds), all design variants show a translation in the direction of the force being applied. The reason for this is the reduced congruence at higher flexion angles, due to the reduced dorsal femoral radii and the reduced compression, as no more body weight is being applied. Likewise, the internal / external rotation follows in the direction of the torque. Here too, the maximum movement can be seen in the swing phase. The Columbus® prosthesis is therefore characterized by high stability in the stance phase and mobility in the swing phase. The motion of the prosthesis allows, among other things, a lower transfer of the forces occurring to the bony implant bed, thereby reducing the risk of early loosening of the implant. The reduced range of motion of Columbus® RP through DD, UC to PS can be explained by the increasing stabilization of the prosthesis designs. Columbus® UC and PS show less mobility in AP and therefore higher stability, enabling the implants to be inserted without PCL.







Figure 27: Kinematic behavior of the Columbus® design variants in the knee wear simulator according to ISO Standard 14243.³²

Columbus® UC



Columbus® PS



Natural kinematics → The motion of the Columbus® prosthesis follows the natural force effect.

Good stability → Deviations of the Columbus® prosthesis are limited.

3. The Columbus® Design

Mobility – Discover stability

3.5 Stability

The high congruence between the femur and gliding surface provides high stability. Laboratory tests show that the Columbus® knee prosthesis withstands stronger stresses than the maximum forces that generally occur in a healthy knee. The risk of subluxation is therefore extremely low with the Columbus® knee prosthesis. The lower stability with regard to torque is explained by the test only being performed in a healthy knee with collateral ligaments. In the test set-up, the function of the collateral ligaments cannot be represented.




Figure 28: Stability of the Columbus® knee prosthesis.³³⁻³⁷

The corresponding gliding surfaces click into the tibial plateaus for the fixed gliding surfaces of Columbus® with a click mechanism and are thereby fixated to the tibial plateau at four fixation points. The PS gliding surfaces are additionally secured with a screw. The fixation of the tibial components has been successfully used for almost three decades, so the risk of subluxation of the gliding surface can be described as extremely unlikely.



Figure 29: Fixation of the fixed variants of the Columbus® tibial plateaus.



Figure 30: Additional securing of the Columbus® PS gliding surfaces with a screw.

Due to the slim wing design, only a small amount of tibial bone needs to be removed for the tibial wing preparation, and the width of the wing provides high rotational stability at the same time.




Figure 31: Design of the Columbus® tibial plateau.

4. Selected materials

Tried and tested – Discover Innovations

4.1 Surface properties

Due to improved surface finishing, the Columbus® knee prostheses achieve an optimum surface quality both for the articulating part of the femoral components and for the PE components.

The abrasion rate can therefore be reduced, particularly in the early post-operative phase as the bedding-in phase of the knee prostheses is significantly reduced.^{38, 39} The surface quality of the fixed tibial plateaus has also been optimized to avoid micro-movements between the tibial plateau and the meniscus components. Wear of the PE on the underside of the meniscus components (backside wear) is therefore minimized.

4.2 Plasmapore® cementless fixation

The cementless implants are available with the established Plasmapore® coating, which has been in widespread use since 1986. The coating is made of pure titanium and is therefore particularly well tolerated. The pore sizes of the Plasmapore® layer are between 50 and 200 µm, which has a positive effect on the growth of the bone tissue.^{42, 43} The very rough Plasmapore® structure primarily supports the fixation of the implants and thereby creates a durable connection between the bone and implant.⁴⁴




Figure 32: Analysis of the surface finish for the Columbus® knee prosthesis according to ISO 7207-2, ASTM F2083.^{40, 41}




Figure 33: Microporous Plasmapore® surface.




Figure 34: Columbus® femoral components with cement pocket and Plasmapore® coating.

4.3 Beta polyethylene

The Columbus® gliding surfaces are manufactured from high-quality Ultra-High-Molecular-Weight Polyethylene (UHMWPE). As defined in ISO 5834-2, the polyethylene is compression molded under low heat. To prevent the risk of oxidation, the implant is packaged in an inert gas atmosphere in the absence of oxygen. Aesculap is convinced of the benefits of moderately cross-linked beta PE and is therefore committed to its use. The gentle beta sterilization procedure considerably improves the oxidation resistance of the polyethylene. The ageing process of the polyethylene is slowed down and the tribological properties of the material are optimized.^{45, 46} The highly cross-linked polyethylene (HXLPE) offered by other prostheses manufacturers show no better abrasion rates than the Columbus® prostheses. There are also no long-term results available for the HXLPE to confirm suitable material properties in knee replacement surgery.^{48, 49}




Figure 35: Oxidation gamma- versus beta-sterilized UHMWPE.

4. Selected materials

Tried and tested – Discover Innovations

4.4 AS Coating

The AS-coated Columbus® prostheses offer an innovative solution for patients who are hypersensitive to metals such as nickel, cobalt and chromium. The multilayer coating is designed to successfully prevent the release of metal ions. The risk of mechanical ablation, which has been known to occur with monolayer coatings, is reduced due to the innovative layer architecture.^{50, 51} The AS coating also improved the wear rate by 60 %.^{52, 53} Almost all the primary Columbus® implant components are available in AS-coated variations.

- **Top layer made from zirconium nitride (approx. 2.5 µm thick)**

The top layer is made from zirconium nitride, a ceramic surface that shows excellent tribological properties. The metal zirconium is one of the five best known metals for biocompatibility.

- **Five intermediate layers made from chromium nitride and chromium carbon nitride (approx. 2.0 µm thick)**

The five intermediate layers are applied alternately. They enable the mechanical integrity of the layer system.

- **Chromium adhesive layer (approx. 0.1 µm thick)**

The adhesive layer of chromium provides a tight connection between the base material and the following layers, thereby creating stability.




Figure 36: Composition of the AS layer architecture.

5. Versatile instruments

Discover precision

5.1 Overview

The Columbus® instrument set offers versatile surgical options.

- Navigated / Manual
- Tibia First / Femur First
- Conventional / MIOS®
- Optional instruments for an individualized surgical procedure.
- "Easy pack" solution for specified surgical procedures.

5.2 OrthoPilot® navigation system

The OrthoPilot® navigation system is proven to help in the precise alignment of the implant with the mechanical axis. The simple and interactive simulation of the femoral and tibial sections, taking into account the individual ligament conditions, enables excellent surgical outcomes.^{54 - 62}



*OrthoPilot®
navigation system*

5. Versatile instruments

Discover precision

In a study of 821 patients, Clemens et al. show that the OrthoPilot® navigation system helps in significantly reducing incorrect positionings of the implant components, thereby increasing the lifetime of the implant.⁶³

Over the last 15 years, more than 150,000 knee prostheses have been successfully implanted using the OrthoPilot® navigation system.




Figure 37: Mechanical axis in 555 navigated vs. 266 manually implanted knee prostheses.⁶³

5.3 Columbus® MIOS® – Minimally Invasive Orthopaedic Solutions

A minimally invasive surgical procedure offers patients the following benefits:

- Smaller incision length
- Less traumatic stress for the soft tissue
- Lower blood loss
- Less bone loss
- Faster post-operative mobilization


Aesculap offers useful instruments for minimally invasive implant surgery:

- Retractors
- Leg holders
- MIOS® saw blades
- MIOS® distractors

Minimally invasive surgical techniques speed up rehabilitation and enable faster recovery of the joint function.^{64, 65} Studies also show that blood loss is reduced when utilizing the minimally invasive surgical technique, and that the hospital stay of the patients can also be reduced.^{66, 67} The narrower approach of the surgical technique, which is kinder on the soft tissue, does however increase what is required of the surgeon. MIOS® instruments were developed with a focus on protecting the soft tissue and they take into account the restricted field of vision. The OrthoPilot® navigation system also enables a more exact positioning of the implant components despite the smaller field of vision.



MIOS® Instruments



Minimally invasive knee replacement surgery with the Columbus® MIOS® instruments.

6. Clinical outcomes

Discover quality

Hauschild et al.⁶⁹

197 CR DD Columbus® knee prostheses were implanted in 125 patients and followed up after 5.8 years.

- The revision rate was 2 %.
- The patients achieved good post-operative stability and kinematics.
- In the 118 cases in which the OrthoPilot® navigation system was used, the implant position was improved.

Picard et al.⁷²

214 cemented Columbus® knee prostheses were implanted in 206 patients with the OrthoPilot® navigation system and followed up after 1 year.

- Problems with anterior knee pain occurred very rarely; in less than 1 % of patients.
- There were only 2 revisions due to infections.
- 98 % of the patients were satisfied or very satisfied.

Walde et al.⁷⁰

60 cementless Columbus® CR knee prostheses were implanted with the OrthoPilot® navigation system and followed up after 5.6 years.

- No prostheses needed to be revised due to aseptic loosening.
- The average KSS increased from 75 to 180 points.
- No patient had problems with anterior knee pain.

Hakki et al.⁷¹

200 cementless Columbus® CR DD knee prostheses were implanted in 182 patients with the OrthoPilot® navigation system and followed up after 2 years.

- The Columbus® knee prosthesis showed good medium-term clinical and radiological outcomes without implant-related problems.
- The patients achieved an average flexion of 122°.
- No prostheses needed to be revised.

Goebel et al.⁷³

109 Columbus® knee prostheses were implanted and compared with 22 NexGen Full Flex knee prostheses (Zimmer, Warsaw, USA) and followed up after an average of 3.9 years.

- In the Knee Society Score Columbus® achieved good results, comparable to those of the NexGen Full Flex prosthesis.
- The operation time was significantly less for Columbus® at 53 minutes (66 minutes for NexGen).
- Columbus® was also judged to be the more bone-sparing implant design and the more flexible knee system.

Malzdorf et al.⁷⁴

116 patients were treated with the Columbus® knee prosthesis. All implant procedures were carried out using the OrthoPilot® navigation system.

- 1 year after the procedure, the patients achieved a flexion angle of 113°.
- Thanks to the OrthoPilot® navigation system, good clinical and radiological outcomes were achieved.

Helwig et al.⁷⁵




193 patients were treated with the Columbus® knee prosthesis and followed up 1 year after the procedure.

- The KSS-F and KSS-K were an average of 170 points one year after the procedure.
- With an OXF of 23 points, high patient satisfaction can be assumed.

Lampe et al.⁷⁶

52 Columbus® prostheses with a fixed platform and 48 Columbus® prostheses with a rotating platform were implanted using the OrthoPilot® navigation system.

- Both the patients with the rotating platform and the patients with the fixed platform showed a good KSS-F 12 months after the procedure with a score of 87 points (before the procedure 47 points) and a good KSS-K with a score of 86 points (before the procedure 30 points).
- The OXF score of 21 points after the procedure (before the procedure 42 points), allows a high level of patient satisfaction to be assumed.



Figures 38: Score results of the Lampe et al. study.

6. Clinical outcomes

Discover quality

Hasegawa et al.⁷⁷

100 Columbus® PS prostheses were implanted using a minimally invasive surgical technique. 50 of these prostheses were implanted using the OrthoPilot® navigation system and 50 prostheses were implanted manually.

- No complications occurred in the post-operative phase.
- The ROM was significantly improved with an average flexion of 121° after 6 months (before the procedure 110°).
- Good post-operative KSS values were achieved both clinically and functionally.
- 94 % of the navigated group fell within the $\pm 3^\circ$ varus / valgus axis deviation. Only 78 % of the conventional group achieved this result ($P = 0.041$).

Stulberg et al.⁷⁸

78 patients were treated with a cementless Columbus® knee prosthesis, 40 of which were manually implanted and 38 with navigation.

- The patients achieved an average flexion of 116.4° six months after the procedure.
- The KSS improved from a pre-operative score of 46.3 points to a score of 84 points 6 months after the procedure.

		MIS Navigated Group	MIS Manual Group	P value
KSS-K	Pre-operative	32 (18)	31 (17)	n.s.*
	Post-operative	95 (7)	92 (9)	n.s.*
KSS-F	Pre-operative	43 (24)	37 (22)	n.s.*
	Post-operative	76 (16)	78 (14)	n.s.*

* n.s.= not significant ($P>0.05$)

Figure 39: Average KSS in the Hasegawa et al. study ($\pm$ standard deviation).





Pre-operative X-ray images



Post-operative X-ray images of the Columbus® CR DD knee prosthesis

7. Important parameters for Columbus® implants

Implant dimensions




Dimensions in (mm)

Size	ML	AP	Box	A	B	C	Pin Z
F1	56	50	34	18.5	34	14	12
F2N	56	53	37	20	36	14.5	13.5
F2	59	53	37	20	36.5	14.5	13.5
F3N	59	56.5	40	21.5	39	16	15
F3	62.5	56.5	40	21.5	39.5	16	15
F4N	62.5	60.5	43.5	23	42	17.5	15
F4	66.5	60.5	43.5	23	42.5	17.5	15
F5N	66.5	65	47.5	26	45.5	20	15
F5	71	65	47.5	26	46	20	15
F6N	71	70	52	28	49	21.5	15
F6	76	70	52	28	49.5	21.5	15
F7	82	75.5	57	30	53.5	23	15
F8	82	80.5	62	32	58	26	15


AP- / ML-dimensions [mm] of the Columbus® femoral implants for any necessary use of intramedullary nails

	AP-nail CR	AP-nail PS	ML-nail
F1	22.5	31	18
F2/F2N	24	32.5	19
F3/F3N	26	34	20.5
F4/F4N	28	36	21
F5/F5N	30	38	22
F6/F6N	32.5	40.5	23
F7	35	42.5	25
F8	39	47	25



Dimensions in (mm)

	T0/T0+	T1/T1+	T2/T2+	T3/T3+	T4/T4+	T5
A	62	65	70	75	80	85
B	41/44	43/46	45/49	48/52	51/55	56
C	14/14.5	15/16	16/17.5	17.5/19	19/20.5	20.5
D	28	28	33	38	43	48
E	12.3	12.3	12.3	12.3	14.3	14.3




Length of the extension stems

Dimensions in (mm)

	T0/T0+	T1/T1+	T2/T2+	T3/T3+	T4/T4+	T5
D	28	28	33	38	43	48
D+S Stem (Short)	80	80	85	90	95	100
D+L Stem (Long)	120	120	125	130	135	140


The overall length of the tibial plateau with the respective extension stem can be calculated by adding Dimension D of the upper table and the stem length Short (52 mm) or Long (92 mm).



The extension stems are available with the diameters 12 and 14 mm.

Patella dimensions

	D _{Patella} x H
Patella P1	Ø 27 mm x 7 mm
Patella P2	Ø 30 mm x 8 mm
Patella P3	Ø 33 mm x 9 mm
Patella P4	Ø 36 mm x 10 mm





Columbus® Implant Matrix – Femoral Parts



Femur CR cemented

Types:	F1	F2N	F2	F3N	F3	F4N	F4	F5N	F5	F6N	F6	F7	F8
Left	NN001K	NN800K	NN002K	NN801K	NN003K	NN899K	NN004K	NN900K	NN005K	NN901K	NN006K	NN007K	NN008K
Left	NN001Z	NN800Z	NN002Z	NN801Z	NN003Z	NN899Z	NN004Z	NN900Z	NN005Z	NN901Z	NN006Z	NN007Z	
Right	NN011K	NN810K	NN012K	NN811K	NN013K	NN909K	NN014K	NN910K	NN015K	NN911K	NN016K	NN017K	NN018K
Right	NN011Z	NN810Z	NN012Z	NN811Z	NN013Z	NN909Z	NN014Z	NN910Z	NN015Z	NN911Z	NN016Z	NN017Z	



Femur PS cemented

Types:	F1	F2N	F2	F3N	F3	F4N	F4	F5N	F5	F6N	F6	F7	F8
Left	NN161K	NN840K	NN162K	NN841K	NN163K	NN939K	NN164K	NN940K	NN165K	NN941K	NN166K	NN167K	NN168K
Left	NN161Z		NN162Z		NN163Z		NN164Z		NN165Z		NN166Z	NN167Z	
Right	NN171K	NN850K	NN172K	NN851K	NN173K	NN949K	NN174K	NN950K	NN175K	NN951K	NN176K	NN177K	NN178K
Right	NN171Z		NN172Z		NN173Z		NN174Z		NN175Z		NN176Z	NN177Z	



Femur CR cementless

Types:	F1	F2N	F2	F3N	F3	F4N	F4	F5N	F5	F6N	F6	F7	F8
Left	NN021K	NN820K	NN022K	NN821K	NN023K	NN919K	NN024K	NN920K	NN025K	NN921K	NN026K	NN027K	NN028K
Right	NN031K	NN830K	NN032K	NN831K	NN033K	NN929K	NN034K	NN930K	NN035K	NN931K	NN036K	NN037K	NN038K

PS Fixation screw for PE in AS version
(generally packed with gliding surface)



Types:	NN497Z	NN498Z	NN499Z
	Height 10/12	Height 14/16	Height 18/20



Patella-3-Peg

Types:	P1	P2	P3	P4
F1-F8	NN481	NN482	NN483	NN484

Columbus® Implant Matrix – Tibial Parts

Tibia cemented											
Types:	T0	T0+	T1	T1+	T2	T2+	T3	T3+	T4	T4+	T5
CR/PS	NN070K	NN058K	NN071K	NN072K	NN073K	NN074K	NN075K	NN076K	NN077K	NN078K	NN079K
CRA/PSA	NN470K	NN469K	NN471K	NN472K	NN473K	NN474K	NN475K	NN476K	NN477K	NN478K	NN479K
RP	–	–	NN271K	NN272K	NN273K	NN274K	NN275K	NN276K	NN277K	NN278K	NN279K
UCR	NN670K	NN668K	NN671K	NN672K	NN673K	NN674K	NN675K	NN676K	NN677K	NN678K	NN679K



CRA/PSA Tibia Hemispacers with screws						
4 mm						
Types:	T0	T1	T2	T3	T4	T5
RM/LL	NN560K	NN563K	NN566K	NN569K	NN572K	NN575K
RL/LM	NN580K	NN583K	NN586K	NN589K	NN592K	NN595K



Tibia cementless											
Types:	T0	T0+	T1	T1+	T2	T2+	T3	T3+	T4	T4+	T5
CR/PS	NN080K	NN059K	NN081K	NN082K	NN083K	NN084K	NN085K	NN086K	NN087K	NN088K	NN089K
RP	–	–	NN281K	NN282K	NN283K	NN284K	NN285K	NN286K	NN287K	NN288K	NN289K

Obturator



Tibia Extension Stems



Types:	T0	T1	T2	T3	T4	T5
RM/LL	NN561K	NN564K	NN567K	NN570K	NN573K	NN576K
RL/LM	NN581K	NN584K	NN587K	NN590K	NN593K	NN596K

Types:	T0	T1	T2	T3	T4	T5
RM/LL	NN561K	NN564K	NN567K	NN570K	NN573K	NN576K
RL/LM	NN581K	NN584K	NN587K	NN590K	NN593K	NN596K

Peek Plug



NN260P

Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN170	NN171	NN172	NN173	–
DD	NN200	NN201	NN202	NN203	NN204	NN205	NN210	NN211	NN212	NN213	NN214
UC fix.	NN400	NN401	NN402	NN403	NN404	NN405	NN410	NN411	NN412	NN413	NN414
UCR	NN600	NN601	NN602	NN603	NN604	NN605	NN610	NN611	NN612	NN613	NN614
RP	–	–	–	–	–	–	NN310	NN311	NN312	NN313	–
PS	NN500	NN501	NN502	NN503	NN504	NN505	NN510	NN511	NN512	NN513	NN514



Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN170	NN171	NN172	NN173	–
DD	NN200	NN201	NN202	NN203	NN204	NN205	NN210	NN211	NN212	NN213	NN214
UC fix.	NN400	NN401	NN402	NN403	NN404	NN405	NN410	NN411	NN412	NN413	NN414
UCR	NN600	NN601	NN602	NN603	NN604	NN605	NN610	NN611	NN612	NN613	NN614
RP	–	–	–	–	–	–	NN310	NN311	NN312	NN313	–
PS	NN500	NN501	NN502	NN503	NN504	NN505	NN510	NN511	NN512	NN513	NN514




Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN140	NN141	NN142	NN143	–
DD	NN130	NN131	NN132	NN133	–	–	NN140	NN141	NN142	NN143	–
UC fix.	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244
UCR	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644
RP	NN330	NN331	NN332	NN333	–	–	NN340	NN341	NN342	NN343	–
PS	NN530	NN531	NN532	NN533	NN534	NN535	NN540	NN541	NN542	NN543	NN544






Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN150	NN151	NN152	NN153	–
DD	NN130	NN131	NN132	NN133	–	–	NN150	NN151	NN152	NN153	–
UC fix.	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244
UCR	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644
RP	NN330	NN331	NN332	NN333	–	–	NN350	NN351	NN352	NN353	–
PS	NN530	NN531	NN532	NN533	NN534	NN535	NN540	NN541	NN542	NN543	NN544









Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN150	NN151	NN152	NN153	–
DD	NN130	NN131	NN132	NN133	–	–	NN150	NN151	NN152	NN153	–
UC fix.	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244
UCR	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644
RP	NN330	NN331	NN332	NN333	–	–	NN350	NN351	NN352	NN353	–
PS	NN530	NN531	NN532	NN533	NN534	NN535	NN540	NN541	NN542	NN543	NN544








Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN150	NN151	NN152	NN153	–
DD	NN130	NN131	NN132	NN133	–	–	NN150	NN151	NN152	NN153	–
UC fix.	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244
UCR	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644
RP	NN330	NN331	NN332	NN333	–	–	NN350	NN351	NN352	NN353	–
PS	NN530	NN531	NN532	NN533	NN534	NN535	NN540	NN541	NN542	NN543	NN544



Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN150	NN151	NN152	NN153	–
DD	NN130	NN131	NN132	NN133	–	–	NN150	NN151	NN152	NN153	–
UC fix.	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244
UCR	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644
RP	NN330	NN331	NN332	NN333	–	–	NN350	NN351	NN352	NN353	–
PS	NN530	NN531	NN532	NN533	NN534	NN535	NN540	NN541	NN542	NN543	NN544



Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN150	NN151	NN152	NN153	–
DD	NN130	NN131	NN132	NN133	–	–	NN150	NN151	NN152	NN153	–
UC fix.	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244
UCR	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644
RP	NN330	NN331	NN332	NN333	–	–	NN350	NN351	NN352	NN353	–
PS	NN530	NN531	NN532	NN533	NN534	NN535	NN540	NN541	NN542	NN543	NN544



Gliding Surfaces											
Types:	T0	T2	T4	T6	T8	T20	T10	T12	T14	T16	T18
CR	–	–	–	–	–	–	NN150	NN151	NN152	NN153	–
DD	NN130	NN131	NN132	NN133	–	–	NN150	NN151	NN152	NN153	–
UC fix.	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244
UCR	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644



Columbus® Implant Matrix – AS Femoral Parts



AS Femur CR cemented

Types:	F1	F2N	F2	F3N	F3	F4N	F4	F5N	F5	F6N	F6	F7	F8
Left	NN001Z	NN800Z	NN002Z	NN801Z	NN003Z	NN899Z	NN004Z	NN900Z	NN005Z	NN901Z	NN006Z	NN007Z	–
Right	NN011Z	NN810Z	NN012Z	NN811Z	NN013Z	NN909Z	NN014Z	NN910Z	NN015Z	NN911Z	NN016Z	NN017Z	–



AS Femur PS cemented

Types:	F1	F2	F3	F4	F5	F6	F7	F8
Left	NN161Z	NN162Z	NN163Z	NN164Z	NN165Z	NN166Z	NN167Z	–
Right	NN171Z	NN172Z	NN173Z	NN174Z	NN175Z	NN176Z	NN177Z	–



Patella-3-Peg

Types:	P1	P2	P3	P4
F1-F8	NN481	NN482	NN483	NN484



Types:	T0	T0+	T1	T1+	T2	T2+	T3	T3+	T4	T4+	T5
CR/PS	NN0702	NN0582	NN0712	NN0722	NN0732	NN0742	NN0752	NN0762	NN0772	NN0782	NN0792
CRA/PSA	NN4702	NN4692	NN4712	NN4722	NN4732	NN4742	NN4752	NN4762	NN4772	NN4782	NN4792
RP	–	–	NN2712	NN2722	NN2732	NN2742	NN2752	NN2762	NN2772	NN2782	NN2792

AS Tibia Extension Stems



short	long	short	long
Ø 12 mm	Ø 12 mm	Ø 14 mm	Ø 14 mm
NN2622	NN2632	NN2652	NN2662

CRA/PSA Tibia Hemispacers with screws



Types:	T0	T1	T2	T3	T4	T5	T0	T1	T2	T3	T4	T5
RM/LL	NN5602	NN5632	NN5662	NN5692	NN5722	NN5752	NN5612	NN5642	NN5672	NN5702	NN5732	NN5762
RL/LM	NN5802	NN5832	NN5862	NN5892	NN5922	NN5952	NN5812	NN5842	NN5872	NN5902	NN5932	NN5962

Obturator



Peek Plug



Ø 12 mm	Ø 14 mm
NN2612	NN2642

Fixation Screw for PS Gliding Surfaces



Types:	NN4972	NN4982	NN4992
	10/12	14/16	18/20

Gliding Surfaces



	T0/T0+												T1/T1+												T2/T2+											
Types:	10	12	14	16	18	20	10	12	14	16	18	20	10	12	14	16	18	20	10	12	14	16	18	20	10	12	14	16	18	20						
CR	—	—	—	—	—	—	NN110	NN111	NN112	NN113	—	—	NN120	NN121	NN122	NN123	—	—	NN120	NN121	NN122	NN123	—	—	NN120	NN121	NN122	NN123	—	—						
DD	NN200	NN201	NN202	NN203	NN204	NN205	NN210	NN211	NN212	NN213	NN214	NN215	NN220	NN221	NN222	NN223	NN224	NN225	NN220	NN221	NN222	NN223	NN224	NN225	NN220	NN221	NN222	NN223	NN224	NN225						
UC fix.	NN400	NN401	NN402	NN403	NN404	NN405	NN410	NN411	NN412	NN413	NN414	NN415	NN420	NN421	NN422	NN423	NN424	NN425	NN420	NN421	NN422	NN423	NN424	NN425	NN420	NN421	NN422	NN423	NN424	NN425						
UCR	NN600	NN601	NN602	NN603	NN604	NN605	NN610	NN611	NN612	NN613	NN614	NN615	NN620	NN621	NN622	NN623	NN624	NN625	NN620	NN621	NN622	NN623	NN624	NN625	NN620	NN621	NN622	NN623	NN624	NN625						
RP	—	—	—	—	—	—	NN310	NN311	NN312	NN313	—	—	NN320	NN321	NN322	NN323	—	—	NN320	NN321	NN322	NN323	—	—	NN320	NN321	NN322	NN323	—	—						
PS	NN500	NN501	NN502	NN503	NN504	NN505	NN510	NN511	NN512	NN513	NN514	NN515	NN520	NN521	NN522	NN523	NN524	NN525	NN520	NN521	NN522	NN523	NN524	NN525	NN520	NN521	NN522	NN523	NN524	NN525						

Gliding Surfaces



	T3/T3+												T4/T4+												T5											
	10	12	14	16	18	20	10	12	14	16	18	20	10	12	14	16	18	20	10	12	14	16	18	20	10	12	14	16	18	20						
CR	NN130	NN131	NN132	NN133	—	—	NN140	NN141	NN142	NN143	—	—	NN150	NN151	NN152	NN153	—	—	NN150	NN151	NN152	NN153	—	—	NN150	NN151	NN152	NN153	—	—						
DD	NN230	NN231	NN232	NN233	NN234	NN235	NN240	NN241	NN242	NN243	NN244	NN245	NN250	NN251	NN252	NN253	NN254	NN255	NN250	NN251	NN252	NN253	NN254	NN255	NN250	NN251	NN252	NN253	NN254	NN255						
UC fix.	NN430	NN431	NN432	NN433	NN434	NN435	NN440	NN441	NN442	NN443	NN444	NN445	NN450	NN451	NN452	NN453	NN454	NN455	NN450	NN451	NN452	NN453	NN454	NN455	NN450	NN451	NN452	NN453	NN454	NN455						
UCR	NN630	NN631	NN632	NN633	NN634	NN635	NN640	NN641	NN642	NN643	NN644	NN645	NN650	NN651	NN652	NN653	NN654	NN644	NN650	NN651	NN652	NN653	NN654	NN644	NN650	NN651	NN652	NN653	NN654	NN644						
RP	NN330	NN331	NN332	NN333	—	—	NN340	NN341	NN342	NN343	—	—	NN350	NN351	NN352	NN353	—	—	NN350	NN351	NN352	NN353	—	—	NN350	NN351	NN352	NN353	—	—						
PS	NN530	NN531	NN532	NN533	NN534	NN535	NN540	NN541	NN542	NN543	NN544	NN545	NN550	NN551	NN552	NN553	NN554	NN555	NN550	NN551	NN552	NN553	NN554	NN555	NN550	NN551	NN552	NN553	NN554	NN555						

Femur stem 6° cemented

Ø 12 mm Ø 15 mm Ø 18 mm



Note: 15 mm femur spacer can be used either distal or posterior. On the other side in the same condyle side, there is no space for

	P0	P1	P2	P3
Types:	Ø 27 mm	Ø 30 mm	Ø 33 mm	Ø 36 mm
FT-E7	NN481	NN482	NN483	NN484

Tibia cemented



with nut inside for ML fixation

Types:	T0	T0+	T1	T1+	T2	T2+	T3	T3+	T4	T4+	T5
	NR070K	NR068K	NR071K	NR072K	NR073K	NR074K	NR075K	NR076K	NR077K	NR078K	NR079K

Tibia offset stems cemented



	Ø 12 mm	Ø 15 mm	Ø 18 mm
Types:	52 mm	92 mm	92 mm
FO-F5	NR191K	NR194K	NR192K
		NR195K	NR193K
			NR196K

Possible Combination Femur / Tibia

Types:	F1	F2	F3	F4	F5	F6	F7
T0/T0+	OK	OK	OK	+	+	+	+
T1/T1+	OK	OK	OK	OK	+	+	+
T2/T2+	–	OK	OK	OK	OK	+	+
T3/T3+	–	–	OK	OK	OK	OK	+
T4/T4+	–	–	+	OK	OK	OK	OK
T5	–	–	+	+	OK	OK	OK

– not compatible
+ possible, but not recommended by Aesculap

Tibia offset stems cementless



	Ø 11 mm	Ø 12 mm	Ø 13 mm	Ø 14 mm	Ø 15 mm
Types:	92 mm	132 mm	92 mm	92 mm	92 mm
T0-T5	NR171K	NR181K	NR172K	NR173K	NR183K
			NR182K	NR178K	NR184K
			NR182K	NR178K	NR185K

	Ø 16 mm	Ø 17 mm	Ø 18 mm	Ø 19 mm	Ø 19 mm
Types:	92 mm	132 mm	92 mm	92 mm	132 mm
T0-T5	NR176K	NR186K	NR177K	NR187K	NR178K
			NR178K	NR188K	NR189K
			NR178K	NR188K	NR190K

PE Gliding Surface Medium Constraint with screw



10 Nm for PE screw

Types:	10 mm	12 mm	14 mm	16 mm	18 mm	20 mm	24 mm	28 mm	32 mm	10 mm	12 mm	14 mm	16 mm	18 mm	20 mm	24 mm	28 mm	32 mm
T0/0+	NR100	NR101	NR102	NR103	NR104	NR105	NR106	–	–	NR600	NR601	NR602	NR603	NR604	NR605	NR606	–	–
T1/1+	NR110	NR111	NR112	NR113	NR114	NR115	NR116	–	–	NR610	NR611	NR612	NR613	NR614	NR615	NR616	–	–
T2/2+	NR120	NR121	NR122	NR123	NR124	NR125	NR126	NR127	–	NR620	NR621	NR622	NR623	NR624	NR625	NR626	NR627	–
T3/3+	NR130	NR131	NR132	NR133	NR134	NR135	NR136	NR137	–	NR630	NR631	NR632	NR633	NR634	NR635	NR636	NR637	–
T4/4+	NR140	NR141	NR142	NR143	NR144	NR145	NR146	NR147	NR148	NR640	NR641	NR642	NR643	NR644	NR645	NR646	NR647	NR648
T5	NR150	NR151	NR152	NR153	NR154	NR155	NR156	NR157	NR158	NR650	NR651	NR652	NR653	NR654	NR655	NR656	NR657	NR658

PE Gliding Surface High Constraint with screw



10 Nm for PE screw

Types:	5 mm	10 mm	15 mm	5 mm	10 mm	15 mm
T0/0+	NR040K	NR041K	NR042K	NR240K	NR241K	NR242K
T1/1+	NR044K	NR045K	NR046K	NR244K	NR245K	NR246K
T2/2+	NR048K	NR049K	NR050K	NR248K	NR249K	NR250K
T3/3+	NR052K	NR053K	NR054K	NR252K	NR253K	NR254K
T4/4+	NR056K	NR057K	NR058K	NR256K	NR257K	NR258K
T5	NR060K	NR061K	NR062K	NR260K	NR261K	NR262K

Tibia Hemi Spacer RM/LL with screw



Tibia Hemi Spacer RL/LM with screw



Literature

- 1 Noble PC, Conditt MA, Cook KF, Mathis KB. The John Insall Award: Patient expectations affect satisfaction with total knee arthroplasty. *Clin Orthop Relat Res.* 2006 Nov;452:35-43.
- 2 Andriacchi TP, Galante JO, Fermier RW. The influence of total knee-replacement design on walking and stair-climbing. *J Bone Joint Surg Am.* 1982 Dec;64(9):1328-35.
- 3 Anouchi YS, McShane M, Kelly F Jr, Elting J, Stiehl J. Range of motion in total knee replacement, *J Clin Orthop Relat Res.* 1996 Oct;(331):87-92.
- 4 Mikashima Y, Tomatsu T, Horikoshi M, Nakatani T, Saito S, Momohara S, Banks SA. In vivo deep-flexion kinematics in patients with posterior-cruciate retaining and anterior-cruciate substituting total knee arthroplasty. *Clin Biomech.* 2010 Jan;25(1):83-7.
- 5 Hakki S, Coleman S, Saleh K, Bilotta VJ, Hakki A. Navigational predictors in determining the necessity for collateral ligament release in total knee replacement. *J Bone Joint Surg Br.* 2009 Sep;91(9):1178-82.
- 6 Delp SL, Kocmond JH, Stern SH. Tradeoffs between motion and stability in posterior substituting knee arthroplasty design. *J Biomech.* 1995 Oct;28(10):1155-66.
- 7 Grupp TM, Schwiesau J. Stability characteristics of the Columbus® UC Total Knee System. Aesculap Test No. V555. Tuttlingen, Germany, 2009 Dec 18.
- 8 Yoshioka Y, Siu D, Cooke TD. The anatomy and functional axes of the femur. *J Bone Joint Surg Am.* 1987 Jul;69(6):873-80.
- 9 Hitt K, Shurman JR 2nd, Greene K, McCarthy J, Moskal J, Hoeman T, Mont MA. Anthropometric measurements of the human knee: correlation to the sizing of current knee arthroplasty systems. *J Bone Joint Surg Am.* 2003;85-A Suppl 4:115-22.
- 10 Mahfouz M, Booth R jr, Argenson J, Merkl BC, Abdel Fatah EE, Kuhn MJ. Analysis of variation of adult femora using sex specific statistical atlases. Paper presented at: 7th International Symposium on Computer Methods in Biomechanics and Biomedical Engineering, Antibes, France, 2006.
- 11 Incavo SJ, Ronchetti PJ, Howe JG, Tranowski JP. Tibial plateau coverage in total knee arthroplasty. *Clin Orthop Relat Res.* 1994 Feb;(299):81-5.
- 12 Mensch JS, Amstutz HC. Knee morphology as a guide to knee replacement. *Clin Orthop Relat Res.* 1975 Oct;(112):231-41.
- 13 Elias SG, Freeman MA, Gokcay EI. A correlative study of the geometry and anatomy of the distal femur. *Clin Orthop Relat Res.* 1990 Nov;(260):98-103.
- 14 Pappas, M.J. and Buechel, F.F. On the use of a constant radius femoral component in meniscal bearing knee replacement, *Journal of Orthopaedic Rheumatology.* 1994;7:27-9.
- 15 Kujala UM, Osterman K, Kormanen M, Nelimarkka O, Hurme M, Taimela S. Patellofemoral relationships in recurrent patellar dislocation. *J Bone Joint Surg Br.* 1989 Nov;71(5):788-92.
- 16 Eckhoff DG, Montgomery WK, Stamm ER, Kilcoyne RF. Location of the femoral sulcus in the osteoarthritic knee. *J Arthroplasty.* 1996 Feb;11(2):163-5.
- 17 Bai B, Baez J, Testa N, Kummer FJ. Effect of posterior cut angle on tibial component loading. *J Arthroplasty.* 2000 Oct;15(7):916-20.
- 18 Postak PD, Heim CS, Greenwald AS. Tibial Plateau Surface Stress in TKA: Factor in Clinical Failure, *Orthopaedic Transactions.* 1993-1994;17(4):1121.
- 19 Postak PD, Heim CS, Greenwald AS. Tibial Plateau Surface Stress in TKA: A factor in Polymer Damage - Series II, *Orthopaedic Transactions.* 1994-1995;18(4):1164.
- 20 McNamara JL, Collier JP, Mayor MB, Jensen RE. A comparison of contact pressures in tibial and patellar total knee components before and after service in vivo. *Clin Orthop Relat Res.* 1994 Feb;(299):104-13.
- 21 Collier JP, McNamara JL, Surprenant VA, Jensen RE, Surprenant HP. All-polyethylene patellar components are not the answer. *Clin Orthop Relat Res.* 1991 Dec;(273):198-203.
- 22 Collier JP, Mayor MB, McNamara JL, Surprenant VA, Jensen RE. Analysis of the failure of 122 polyethylene inserts from uncemented tibial knee components. *Clin Orthop Relat Res.* 1991 Dec;(273):232-42.
- 23 Bartel DL, Bicknell VL, Wright TM. The effect of conformity, thickness, and material on stresses in ultra-high molecular weight components for total joint replacement. *J Bone Joint Surg Am.* 1986 Sep;68(7):1041-51.

- 24 Robertsson O, Knutson K, Lewold S, Lidgren L. The Swedish Knee Arthroplasty Register 1975-1997: an update with special emphasis on 41,223 knees operated on in 1988-1997. *Acta Orthop Scand*. 2001 Oct;72(5):503-13.
- 25 Willert HG, Bertram H, Buchhorn GH. Osteolysis in alloarthroplasty of the hip. The role of ultra-high molecular weight polyethylene wear particles. *Clin Orthop Relat Res*. 1990 Sep;(258):95-107.
- 26 Schwiesau J. Determination of the wear behaviour of the Columbus® PS Knee System. Aesculap Test No. T012. Tuttlingen, Germany, 2007 Mar 08.
- 27 Grupp TM, Kaddick C, Schwiesau J, Maas A, Stulberg SD. Fixed and mobile bearing total knee arthroplasty-influence on wear generation, corresponding wear areas, knee kinematics and particle composition. *Clin Biomech (Bristol, Avon)*. 2009 Feb; 24(2):210-7. 2008 Dec 31.
- 28 Kretzer JP, Jakubowitz E, Hofmann K, Heisel C, Lietz E. Influence of resting periods on polyethylene wear in a knee simulator study. *Med Sci Monit*. 2009 Nov;15(11):MT143-6.
- 29 Kretzer JP, Jakubowitz E, Sonntag R, Hofmann K, Heisel C, Thomsen M. Effect of joint laxity on polyethylene wear in total knee replacement. *J Biomech*. 2010 Apr 19;43(6):1092-6.
- 30 Schwiesau J. Determination of the wear behaviour of the Columbus® UCR Knee System. Aesculap Test No. T63. Tuttlingen, Germany, 2010 Sep 09.
- 31 Grupp TM, Stulberg D, Kaddick C, Maas A, Fritz B, Schwiesau J, Bloemer W. Fixed bearing knee congruency - influence on contact mechanics, abrasive wear and kinematics. *Int J Artif Organs*. 2009 Apr;32(4):213-23.
- 32 Schwiesau J, Grupp TM, Hintner M, Kaddick C, Hermle T, Blomer W. Auswirkungen unterschiedlicher Kopplungsgrade einer Knieprothese auf Kinematik und Verschleiss, Ergebnisse einer Simulatorstudie nach ISO 14243-1,2, Internationale Biomechanik- und Biomaterial-Tage München, 2004.
- 33 Seireg A, Arvikar RJ. The Prediction of the muscular load sharing and joint Forces in the lower extremities during Walking, *J of Biomech*. 1975;8:89-102.
- 34 Postak PD, Greenwald AS. Stability Characteristics of the Columbus® CR Knee System. Aesculap Test No. V449. T023J. Tuttlingen, Germany, 2002.
- 35 Postak PD, Greenwald AS. Stability Characteristics of the Columbus® Deep Dish Knee System. Aesculap Test No. V450. Tuttlingen, Germany, 2002.
- 36 Postak PD, Greenwald AS. Stability Characteristics of the Columbus® PS Knee System. Aesculap Test No. V451. Tuttlingen, Germany, 2003.
- 37 Postak PD, Greenwald AS. Stability Characteristics of the Columbus® UC Knee System. Aesculap Test No. V555. Tuttlingen, Germany, 2006.
- 38 Ingham E, Fisher J. Biological reactions to wear debris in total joint replacement. *Proc Inst Mech Eng H*. 2000;214(1):21-37.
- 39 Que L, Topoleski LD, Parks NL. Surface roughness of retrieved CoCrMo alloy femoral components from PCA artificial total knee joints. *J Biomed Mater Res*. 2000;53(1):111-8.
- 40 Grupp TM, Schneid U. Surface Roughness Analysis of the Columbus® Knee System according to ISO 7207-2 and ASTM F2083-06. Aesculap Test No. V557, Tuttlingen, Germany, 2007.
- 41 Grupp TM, Holeried M. Surface Roughness Analysis of the Columbus® UCR Tibia Plateau component NN675K and the Columbus® UCR Bearing component NN630 according to ISO 7207-2:1998 (E) and ISO 21536:2007 (E). Aesculap Test No. V728. Tuttlingen, Germany, 2010 June.
- 42 Predecki P, Stephan JE, Auslaender BA, Mooney VL, Kirkland K. Kinetics of bone growth into cylindrical channels in aluminum oxide and titanium. *J Biomed Mater Res*. 1972 Sep;6(5):375-400.
- 43 Bobyn JD, Pilliar RM, Cameron HU, Weatherly GC. The optimum pore size for the fixation of porous-surfaced metal implants by the ingrowth of bone. *Clin Orthop Relat Res*. 1980 Jul-Aug;(150):263-70.
- 44 Pitto RP, Böhner J, Hofmeister V. Factors affecting the primary stability of acetabular components. An in vitro study. *Biomed Tech*. 1997 Dec;42(12):363-8.
- 45 Bell CJ, Walker PS, Abeysondera MR, Simmons JM, King PM, Blunn GW. Effect of oxidation on delamination of ultrahigh-molecular-weight polyethylene tibial components. *J Arthroplasty*. 1998 Apr;13(3):280-90.

Literature

- 46 Blömer W, Lohrmann E. Verschleißbeständigkeit von UHMWPE-Artikulationen in der Hüftgelenkendoprothetik. In: Weller S, Braun A, Eingartner C, Maurer F, Weise K, Winter E, Volkmann R, editors. Das Bicontact® Hüftendoprothesensystem 1987-2007. Stuttgart: Georg Thieme Verlag; 2007:94-100.
- 47 Kurtz SM, Muratoglu OK, Evans M, Edidin AA. Advances in the processing, sterilization, and crosslinking of ultra-high molecular weight polyethylene for total joint arthroplasty. *Biomaterials*. 1999 Sep;20(18):1659-88.
- 48 Muratoglu OK, Mark A, Vittetoe DA, Harris WH, Rubash HE. Polyethylene damage in total knees and use of highly crosslinked polyethylene. *J Bone Joint Surg Am*. 2003;85-A Suppl 1:S7-S13.
- 49 Sychterz CJ, Shah N, Engh CA. Examination of wear in Duraloc acetabular components: two- to five-year evaluation of Hylamer and Enduron liners. *J Arthroplasty*. 1998 Aug;13(5):508-14.
- 50 Reich J, Hovy L, Lindenmaier HL, Zeller R, Schwiesau J, Thomas P, Grupp TM. Präklinische Ergebnisse beschichteter Knieimplantate für Allergiker. *Orthopade*. 2010 May;39(5):495-502.
- 51 Thull R, Gradinger R. Allergy Solution for orthopaedic knee implants. *Biomaterialien* 2007;8(2):109-10.
- 52 Grupp TM, Affatato S, Jansson V, Müller P, Schwiesau J, Fritz B, Blömer W. Biotribology of a mobile unicompartimental knee bearing design. 18th European Orthopaedic Research Society Meeting, Davos, Switzerland 2010.
- 53 Affatato S, Spinelli M, Lopomo N, Grupp TM, Marcacci M, Toni A. Can the method of fixation influence the wear behaviour of ZrN coated unicompartimental mobile knee prostheses? *Clin Biomech* 2011 Feb;26(2):152-8.
- 54 Decking R, Markmann Y, Fuchs J, Puhl W, Scharf HP. *J Arthroplasty*. 2005 Apr;20(3):282-8. Leg axis after computer-navigated total knee arthroplasty: a prospective randomized trial comparing computer-navigated and manual implantation. *J Arthroplasty*. 2005 Apr;20(3):282-8.
- 55 Hart R, Janecek M, Chaker A, Bucek P. Total knee arthroplasty implanted with and without kinematic navigation. *Int Orthop*. 2003;27(6):366-9.
- 56 Jenny JY, Boeri C. Computer-assisted implantation of a total knee arthroplasty: a case-controlled study in comparison with classical instrumentation. *Rev Chir Orthop Reparatrice Appar Mot*. 2001 Nov;87(7):645-52.
- 57 Jenny JY, Clemens U, Kohler S, Kiefer H, Konermann W, Miehke RK. Consistency of implantation of a total knee arthroplasty with a non-image-based navigation system: a case-control study of 235 cases compared with 235 conventionally implanted prostheses. *J Arthroplasty*. 2005 Oct;20(7):832-9.
- 58 Saragaglia D, Picard F, Chaussard C, Montbarbon E, Leitner F, Cinquin P. Computer-assisted knee arthroplasty: comparison with a conventional procedure. Results of 50 cases in a prospective randomized study. *Rev Chir Orthop Reparatrice Appar Mot*. 2001 Feb 1;87(1):18-28.
- 59 Haaker RG, Stockheim M, Kamp M, Proff G, Breitenfelder J, Ottersbach A. Computer-assisted navigation increases precision of component placement in total knee arthroplasty. *Clin Orthop Relat Res*. 2005 Apr;(433):152-9.
- 60 Song EK, Seon JK, Yoon TR, Park SJ, Bae BH, Cho SG. Functional results of navigated minimally invasive and conventional total knee arthroplasty: a comparison in bilateral cases. *Orthopedics*. 2006 Oct;29(10 Suppl):145-7.
- 61 Mason JB, Fehring TK, Estok R, Banel D, Fehring K. Meta-analysis of alignment outcomes in computer-assisted total knee arthroplasty surgery. *J Arthroplasty*. 2007 Dec;22(8):1097-106.
- 62 Jeffery RS, Morris RW, Denham RA. Coronal alignment after total knee replacement. *J Bone Joint Surg Br*. 1991 Sep; 73(5):709-14.
- 63 Clemens U, Miehke RK, Kohler S, Kiefer H, Jenny JY, Konermann W. Computer assisted navigation with the Ortho-Pilot System and the SEARCH Evolution prosthesis - Results of a multicenter Study. In: Konermann W, Haaker R, editors. *Navigation und Robotic in der Gelenk- und Wirbelsäulenchirurgie*. Heidelberg: Springer-Verlag; 2003:207-16.
- 64 Bohlen KE, Dries SPM, Sufi-Siavach A, Hille E, Lampe F. Quality of implant alignment and results in minimally invasive navigated total knee arthroplasty. 6th annual conference of the international Society for Computer Assisted Orthopaedic Surgery, Montréal, Quebec Canada. June 21-24, 2006.
- 65 Song EK, Seon JK, Yoon TR, Park SJ, Bae BH, Cho SG. Functional results of navigated minimally invasive and conventional total knee arthroplasty: a comparison in bilateral cases. *Orthopedics*. 2006 Oct;29(10 Suppl):S145-7.

- 66 Hart R, Janecek M, Cizmár I, Stipčák V, Kucera B, Filan P. Minimally invasive and navigated implantation for total knee arthroplasty: X-ray analysis and early clinical results. *Orthopade*. 2006 May;35(5):552-7.
- 67 Coon TM. The economic impact of minimally invasive total knee arthroplasty. *Am J Orthop*. 2006 Jul;35(7 Suppl):33-5.
- 68 Grupp TM, Pfaff A, Miranda MK. Columbus® MIOS® Tibia Tray (CR, DD, UC, PS) primary stability in natural bone; comparison to Columbus® Standard Tray. Aesculap Test No. I010. Tuttlingen, Germany, Nov 15, 2006.
- 69 Hauschild O, Thielemann F. Mid-term results of the Columbus® knee endoprosthesis system. Submitted to: SICOT Meeting 2011. Prague, Sep 6-9 2011.
- 70 Schüttrumpf JP, Balcerek P, Wachowski MM, Stürmer KM, Walde H-J, Walde TA. Navigated Cementless Total Knee Arthroplasty – Medium-Term Clinical and Radiological Results. Submitted to: *Open Orthopedics Journal*. Mar 2012. Walde TA, Burgdorf D, Walde HJ. Process optimization in navigated total knee arthroplasty. *Orthopedics*. 2005 Oct;28 (10 Suppl):1255-8.
- 71 Hakki S. A five year experience with OrthoPilot® navigated Columbus® total knee replacement. Paper presented at: CHAOS, 8th Annual meeting of the international society for Computer Assisted Surgery, Hong Kong, 2008.
- 72 Katipalli G, Deakin AH, Greaves J, Reynolds T, Picard F. One-year followup of 214 total knee arthroplasties with navigated columbus implants. *Orthopedics*. 2008 Oct;31(10 Suppl 1).
- 73 Goebel D, Schultz W. The Columbus® Knee system: A bone-preserving alternative to high flexion implants? *Archives of Orthopaedic and Trauma Surgery*. Dec 2011; Ausgabe: 9-5.
- 74 Becker R, Malzdorf M, Stärke C, Randolph P, Lohmann C. No difference between tibia-first and femur-first techniques in TKA using computer-assisted surgery. *Knee Surg Sports Traumatol Arthrosc*. 2012 Feb 25.
- 75 Helwig P, Ateschrang A, Hömig D, Müller P, Knäse W. Multicenterstudie Columbus®. Unpublished Aesculap Internal Study Report, 2011.
- 76 Lampe F, Bohlen KM, Dries SPM, Sufi-Siavach A, Hille E. One year after navigated total knee replacement, no clinically relevant difference between fixed bearing and mobile bearing knee replacement in a double-blind randomized controlled trial. *Open orthop J*. 2011;5:201-8. 2011 May.
- 77 Hasegawa M, Yoshida K, Wakabayashi H, Sudo A. Minimally invasive total knee arthroplasty: comparison of jig-based technique versus computer navigation for clinical and alignment outcome. *Knee Surg Sports Traumatol Arthrosc*. 2011 Jun;19(6):904-10.
- 78 Stulberg SD, Yaffe MA, Koo SS. Computer assisted Surgery vs. manual Total knee arthroplasty: a case controlled study. *J Bone Joint Surg*. 2006;88A(4):47-54.

